


Manufactured by MACHEREY-NAGEL  
**OPTIMA<sup>®</sup>**  
**LIPODEX<sup>®</sup>**  
**HYDRODEX**  
Made in Germany



User Manual



# GCC

Fused silica capillary column

**MACHEREY-NAGEL**



## Fused Silica Kapillarsäulen für die GC · Einbau- und Gebrauchsanleitung

Installation der Säule

1

Inbetriebnahme und Test der Säule und des Chromatographen

3

Allgemeine Informationen zur Säulenhandhabung

6

Hilfe

7

## Fused silica capillary columns for GC · Installation and operating instructions

Column installation

8

Start-up and testing column and chromatograph

10

General information for handling

13

Technical support and customer services

14

Visit our chromatography pages

[www.mn-net.com/chroma](http://www.mn-net.com/chroma)

- Free application database with more than 3000 chromatography applications
- Derivatization protocols
- Integrated webshop
- Product information
- Safety data sheets
- ... and much more

### Please read these instructions BEFORE installing – and certainly before heating – your column!

The present capillary column from MACHEREY-NAGEL has been manufactured in our own facilities and will be the heart of your gas chromatographic analytical system. It has been produced following state-of-the-art technologies and has been individually tested; the enclosed certificate documents its quality standard.

In order to obtain good results from your column as long as possible, you should carefully read and follow these instructions. We have grouped them into 4 chapters:

- I Column installation
- II Start-up and testing column and chromatograph
- III General information for handling
- IV Technical support and customer services

### I. Column installation

The following items are required for proper column installation:

- Ferrules (with hole appropriate for the respective column outer diameter)

Column ID [mm]	Column OD [mm] (= bore of ferrule)
0.05–0.20	0.40
0.25	0.40
0.32	0.50
0.53	0.80

For GC instruments of Agilent special ferrules are available – please see our website [www.mn-net.com](http://www.mn-net.com) or our catalog Chromatography.

- Diamond file (e.g., REF 708300)
- Magnifying lens (e.g., REF 706296)
- Vial with volatile solvent (e.g., pentane, cyclohexane, acetone)
- Leak detector or leak check liquid
- Flow meter

**Important note:** Do not remove the cords from the capillary column! The cords attach the column to the cage and prevent unnecessary contact between the column and metal parts of the cage. The cords are temperature stable.

**Caution:** Without carrier gas flow at high temperatures any stationary phase in a capillary column will be destroyed in no time!


Now you are ready for installation:

**Caution:** Avoid sharp bends. The column should be installed in an unstrained configuration. A column hanger should only contact the cage. For coupling of columns use only zero dead volume fittings or devices to avoid peak broadening, tailing or even memory effects.

- a) Open the plastic packing (cut – not tear) and remove the column.
- b) Check the column for mechanical transport damage; there must be no free ends visible (ends not melted)! If this should be the case, return the column to your distributor for replacement.
- c) The melted column ends have to be removed. For this purpose, score the polyimide layer about 5 cm from the melted end using a diamond file and break the end off straight. Make sure that no column fragments get into the column. For this purpose it is advisable to take off the removed end in a downward movement.
- d) Unwind both column ends from the cage for about one turn.
- e) Suspend the column in the middle of the oven.
- f) Place a column nut and a ferrule over each column end, then cut about 5 cm from each column end using a diamond file (in order to remove any ferrule particles from the column).

## Fused silica capillary column for GC

- g) Examine the column ends with a magnifying lens (REF 706296). If the column ends are not smooth or even show strong burrs, cut off another piece from the column.



- h) Consult the manual of your GC instrumentation because it is very important to install the column entrance at the correct insertion distance in the injector. Insert one column end into the injector and attach it according to the instructions of the instrument manufacturer, i.e., fix it only slightly more than finger-tight.
- i) Adjust the flow of the carrier gas through the injector according to the instrument manual, the analytical instructions or the average values of the following table:

ID	Length	Pressure	Approximate dead time
[mm]	[m]	[bar]	[min]
0.05	10	2.5 (H <sub>2</sub> )	1.5
	20	4.5 (H <sub>2</sub> )	3
0.20	12	0.8 (He)	1.0
	25	1.5 (He)	1.8
	50	3.0 (H <sub>2</sub> )	3.0
0.25	10	0.3 (He)	1.0
	15	0.4 (He)	1.2
	25	0.7 (He)	2.5
	30	0.8 (He)	2.8
	50	1.5 (He)	5
	60	1.8 (He)	5.5
0.32	10	0.3 (He)	0.8
	15	0.4 (He)	0.9
	25	0.5 (He)	1.4
	30	0.6 (He)	1.6
	50	1.0 (He)	3.5
	60	1.2 (He)	4.0
0.53	10	0.2 (He)	0.6
	15	0.2 (He)	0.7
	25	0.3 (He)	1.0
	30	0.4 (He)	1.1
	50	1.0 (H <sub>2</sub> )	1.5

- j) Place the free end of the capillary column into a vial with a volatile solvent. After not more than 30 s gas bubbles should indicate a steady flow through the column. If no bubbles appear, repeat steps f) – i). If you still don't see a gas flow check the following:

- the supply of carrier gas
  - the pressure of the carrier gas in the supply pipes (pressure regulator, valves)
  - if the capillary column is broken
- k) Insert the second column end into the detector and attach it according to the instructions of the instrument manufacturer, i.e., fix the nut slightly. Exception: If using highly sensitive detectors (MS, ECD) you should first condition the column (see below) and then connect it to the detector.

*Hint: For subsequent column installations you may mark the installation depth of the column in the injector and detector in the gas chromatograph using a moisture-proof felt tip pen (e.g., on the door or the edge of the oven).*

- l) Check the column connections in the oven for leaks (e.g., with leak check liquid; this is especially important when using hydrogen as carrier gas)!
- m) Close the oven door, start the gas chromatograph:
- Heat the injector and the detector (maximum to the long-term temperature limit of the capillary column) and ignite, if required (e.g., FID).
  - Heat the oven (choose a temperature between 60 and 100 °C).
  - Start the integrator or the data processing system.
  - Adjust the split ratio (typical 1:30).
  - Inject methane or a volatile solvent to check the overall performance of the gas-chromatographic system.
- n) Before first usage and as required, condition the column, (see chap. II.1).
- o) Possibly reproduce the column test from certificate (see chap. II.2) – now the column is ready for use.

## II. Start-up and testing column and chromatograph

### 1. Conditioning

Each column has been conditioned prior to its quality control. Conditioning is used to elute volatile components of sample molecules and stationary phase and to remove active sites of the inner column surface.

Conditioning should always be performed after a longer storage of the column and performed until the baseline is stable (usually approx. 4 h).

The two numerical temperature values mentioned in the certificate indicate the (lower) isothermal long-term temperature (= conditioning temperature) and the (higher) short-term temperature (= maximum in temperature gradients during 5–15 min).

During conditioning never exceed the long-term temperature, and make sure that the gas supply is never interrupted.

### 2. Column test

Each column from MACHEREY-NAGEL is tested in a gas chromatograph, and the individual chromatogram with conditions (the test certificate) is supplied with the column. It should be retained for future reference.

If possible you should reproduce the included test chromatogram before using the column for your analytical tasks. The respective test mixture can be separately purchased (see [www.mn-net.com](http://www.mn-net.com)).

This will indicate whether

- column installation is correct
- samples are properly injected
- the overall chromatographic system including evaluation works accurately.

By reproducing the column test you document the performance of the column, and thus you can at any time compare the state of the column with the initial quality.

Using the column test, you can identify a number of typical failures; this is the reason why reproduction of this test should always be the first step in troubleshooting.

If the test result is satisfactory, the planned application can be operated.

### 3. Checking column installation

What should you do when your own test does not show the same result as the test certificate?


- a) Check the correct installation of the column and test again.
- b) Compare the faulty chromatogram with the examples in our troubleshooting and try to correct the problem.
- c) Contact your nearest MACHEREY-NAGEL office or your local distributor and ask for help.

### 4. Troubleshooting




If you cannot reproduce the test chromatogram of your capillary column, this may be due to one of the following reasons (see page 11 and 12).

For example, a dead volume in the injector produces tailing. This is especially strong for early eluting components in a temperature program. Dead volume in the detector also results in tailing, however, in this case all components are affected.

## Fused silica capillary column for GC

<i>Observation</i>	<i>Possible cause</i>	<i>Suggested remedy</i>
<p><b>Shortcoming A</b></p>  <p><b>Missing or overlapping peaks</b></p>	<p>Sample too diluted</p> <p>Column connection leaks, column not properly installed</p> <p>Injector temperature too low</p> <p>Sample decomposes in the injector</p> <p>Column adsorbs substances or decomposes them</p>	<p>Increase injection volume, concentrate sample</p> <p>Check for leaks, change ferrules, reinstall column</p> <p>Check temperature program, increase injector temperature</p> <p>Check temperature program, reduce injector temperature, change liner, check capillary ends</p> <p>Check capillary ends, check intact deactivation using a suitable test mixture; for poor results shorten both column ends by about 10 cm; or replace column</p> <p>If the column test does not show any defects: a) use a column with thicker film; b) use alternative phase</p>
<p><b>Shortcoming B</b></p>  <p><b>Increasing baseline, at high temperatures bleeding, strong noise</b></p>	<p>Septum particles in the column entrance</p> <p>Heating program too fast, final temperature too high</p> <p>Column bleeding due to incomplete column conditioning</p> <p>Soiled detector</p>	<p>Cut off one turn from the column entrance, change injection septa</p> <p>Decrease temperature gradient and final temperature</p> <p>Disconnect detector and condition column (also see instructions for your GC instrumentation and detector)</p> <p>Clean detector according to manufacturers instructions</p>
<p><b>Shortcoming C</b></p>  <p><b>Regular interfering peaks</b></p>	<p>Bleeding of silicone septa</p>	<p>Use new septum or septum grade with lower bleeding (high temperature septum)</p>
<p><b>Shortcoming D</b></p>  <p><b>Irregular interfering peaks / spikes / ghost peaks</b></p>	<p>Impurities from sample vials and/or sample preparation</p> <p>Dirty syringe</p> <p>Sample is decomposed</p> <p>Column adsorbs substances or decomposes them</p> <p>Bleeding of butyl septa</p>	<p>Check sample preparation and autosampler vials, use septa with lower bleeding</p> <p>Use clean syringe</p> <p>Check temperature program, oven temperature (thermometer); when chromatographing thermolabile compounds, reduce injector temperature, change liner</p> <p>Check capillary ends, check intact deactivation using a suitable test mixture; for poor results shorten both column ends by about 10 cm; or replace column</p> <p>If the column test does not show any defects: a) use a column with thicker film; b) use alternative phase</p> <p>Use silicone septa for injector and sample vials</p>

## Installation and operating instructions

<i>Observation</i>	<i>Possible cause</i>	<i>Suggested remedy</i>
<b>Shortcoming E</b> 	Column overload Sample is decomposed  Column adsorbs substances or decomposes them	Reduce injection volume Check temperature program, oven temperature (thermometer); when chromatographing thermolabile compounds, reduce injector temperature, change liner Check capillary ends, check intact deactivation using a suitable test mixture; for poor results shorten both column ends by about 10 cm; or replace column  If the column test does not show any defects: a) use a column with thicker film; b) use alternative phase
<b>Shortcoming F</b> 	Sample with high boiling point System leaks Sample is decomposed  Column adsorbs substances or decomposes them	Derivatize polar, basic or high-boiling samples Check for proper column installation and leaks, change ferrules Check temperature program, oven temperature (thermometer); when chromatographing thermolabile compounds, reduce injector temperature, change liner Check capillary ends, check intact deactivation using a suitable test mixture; for poor results shorten both column ends by about 10 cm; or replace column  If the column test does not show any defects: a) use a column with thicker film; b) use alternative phase
<b>Tailing, heavy peak broadening of the peak back</b>	Compounds which always show strong tailing	No remedy possible
<b>Shortcoming G</b> 	Split flow too low Column overload	Increase split flow Inject less or increase split flow
<b>Broad peaks</b>		

### III. General information for handling

#### 1. Column regeneration

If contamination is the cause of a loss in column performance (compare shortcoming D) you may try one of the following remedies:

- a) Remove the contaminated end from the column, i.e. in general shorten the injector column end by about 50 cm (at least one column turn).
- b) Remove column end in the transfer line of the detector, because here phase deterioration may be faster than in the rest of the column due to the permanent high temperatures.
- c) You may also try to "flood" the column in the GC using 20 µL acetone followed by a normal heating program.
- d) The column can be rinsed (*only chemically bonded phases*). For this purpose disconnect the capillary column from the gas chromatograph and connect the detector end (mark it!) to the rinsing device. For rinsing use 10 mL pentane, acetone or toluene for nonpolar – medium polar phases or dichloromethane for polar phases. After rinsing dry the stationary phase for at least 2 h in a stream of carrier gas at ambient temperature before you install the column in the GC and heat it in a temperature program (1 °C/min) up to the maximum isothermal temperature. Hold this temperature for about 1 h.

#### 2. Column storage

If you want the capillary column to stay in the GC over night we recommend to reduce the split flow and adjust the oven temperature to 80–100 °C.

If you have to turn off the carrier gas or if it is necessary to change the column, the oven should be cooled to ambient temperature before disconnecting the column. In order to prevent oxygen damaging of the stationary phase the column should be stored with sealed column ends (*especially important for Wax and FFAP columns*). Fused silica columns can be easily sealed by sticking the ends (sideways) into a septum. For long-term storage we recommend flame sealing and storage in the original shipping container.

#### 3. Maintenance of value

The column at hand consists of a tube of fused silica manufactured with high precision and coated with a protective polyimide layer which results in a high column flexibility.

**Caution:** It is very important to avoid damage to the polyimide layer, since unprotected fused silica becomes brittle and thus easily breaks.

The column is coiled onto a cage which largely protects it from mechanical damage. Avoiding mechanical or extreme thermal strain (e.g., opening the oven door at high temperatures!).

In addition, you can protect your column from chemical stress by performing a suitable sample preparation, by controlling the pH value of your sample (neither acid nor basic) and by using a precolumn. If matrix problems occur, a change in sample preparation can considerably extend the lifetime of the column.

We recommend that you use CHROMAFIL® syringe filters for sample clarification, CHROMABOND® SPE columns for selective sample preparation and/or deactivated capillary columns (precolumns) to protect the separation column from contamination. If necessary 0.5–1 m (always at least one column turn) of the guard column can be removed without need to shorten the analytical column.



### IV. Technical support and customer services

If you have any questions concerning the use of capillary columns please contact us:

MACHEREY-NAGEL GmbH & Co. KG

Germany and international	Phone	+49 24 21 969-0
	Fax	+49 24 21 969-199
	E-mail	info@mn-net.com
Switzerland	Phone	+41 62 388 55 00
	Fax	+41 62 388 55 05
	E-mail	sales-ch@mn-net.com
France	Phone	+33 388 68 22 68
	Fax	+33 388 51 76 88
	E-mail	sales-fr@mn-net.com
USA	Phone	+1 484 821 0984
	Fax	+1 484 821 1272
	E-mail	sales-us@mn-net.com

### Liability

MACHEREY-NAGEL products are intended for general laboratory use only. They are not suited for clinical use (diagnostic, therapeutic or prognostic). MACHEREY-NAGEL does not assume any responsibility for damages due to improper application of our products or due to application of our products in other fields of application. The user has to ensure that the products used are suitable for the planned application. Application on the human body is strictly forbidden! The user is liable of all damages resulting from such application.

MACHEREY-NAGEL products are suited for qualified personnel only. All goods have to be checked immediately upon receipt. Damage claims are only acceptable in writing within 8 working days of receipt of all the goods. In case of legitimate claims the buyer can only require replacement of the goods. If replacement is not possible, the buyer has the right to choose between alternative products with same value or refund. The buyer cannot claim further compensation. All returns must first be authorized by us in writing. We offer a one-year warranty from date of delivery that our products will conform to applicable specifications set forth in the product specifications if not sold to persons set forth in § 13 of German Civil Code. In such a case, the provisions of the German Civil Code shall be valid.

We will be glad to help you with:

- Technical support by phone
- Technical documentation as comprehensive collections of applications for our products
- Local technical support (e.g., by your local distributor)
- Our extensive internet collection of applications for GC, HPLC, TLC, SPE
- Visit our chromatography pages: [www.mn-net.com](http://www.mn-net.com)


This warranty does not cover defects in or damage to the products which are due to improper installation or maintenance, misuse, neglect or any use other than ordinary commercial application.

Any discharge from liability will be void if a defect results from a heavily negligent or intentional breach of contract or if the buyer will be bodily injured as a result from a negligent or intentional breach.

The seller shall not be liable for the goods being fit for a particular purpose unless otherwise agreed upon, to which the buyer intends to put them. This warranty is strictly exclusive. Any further damage compensation is impossible.

MACHEREY-NAGEL makes no other warranty of any kind whatsoever, and SPECIFICALLY DISCLAIMS AND EXCLUDES ALL OTHER WARRANTIES OF ANY KIND OR NATURE WHATSOEVER; DIRECTLY OR INDIRECTLY; EXPRESS OR IMPLIED; INCLUDING; WITHOUT LIMITATION; AS TO THE SUITABILITY; REPRODUCIBILITY; DURABILITY; FITNESS FOR A PARTICULAR PURPOSE OR USE; MERCHANTABILITY; CONDITION; OR ANY OTHER MATTER WITH RESPECT TO MACHEREY-NAGEL PRODUCTS OR THE SALE OF MACHEREY-NAGEL PRODUCTS.

If not superseded by these terms the terms and conditions for the sale of goods (export version) are additionally valid.



... we Meet your Needs

**Quality "Made in Germany"  
For research purposes only**

Experience our professional  
technical support in chromatography

E-mail: [info@mn-net.com](mailto:info@mn-net.com)

Tel.: + 49 24 21 969-0

Visit our online application database at

[www.mn-net.com/apps](http://www.mn-net.com/apps)



GC manual Rev08 2/5/07.16 PD  
A033487  
Printed in Germany



MACHERY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6–8 · 52355 Düren · Germany

DE / International:

Tel.: +49 24 21 969-0

Fax: +49 24 21 969-199

E-mail: [info@mn-net.com](mailto:info@mn-net.com)

CH:

Tel.: +41 62 388 55 00

Fax: +41 62 388 55 05

E-mail: [sales-ch@mn-net.com](mailto:sales-ch@mn-net.com)

FR:

Tel.: +33 388 68 22 68

Fax: +33 388 51 76 88

E-mail: [sales-fr@mn-net.com](mailto:sales-fr@mn-net.com)

US:

Tel.: +1 484 821 0984

Fax: +1 484 821 1272

E-mail: [sales-us@mn-net.com](mailto:sales-us@mn-net.com)