

# TABELA DE ABERTURAS DE TELAS PARA PENEIRAS

ISO / ASTM	MESH / TYLER	Abert. mm.	TOL. ( malha)	
			DE	ATÉ
		<b>Nominal</b>		
4"	****	<b>100,00</b>	<b>97,06</b>	<b>102,94</b>
3.1/2"	****	<b>90,00</b>	<b>87,34</b>	<b>92,66</b>
3"	****	<b>75,00</b>	<b>72,78</b>	<b>77,22</b>
2.1/2"	****	<b>63,00</b>	<b>61,13</b>	<b>64,87</b>
2"	****	<b>50,00</b>	<b>48,51</b>	<b>51,49</b>
1.3/4"	****	<b>45,00</b>	<b>43,65</b>	<b>46,35</b>
1.1/2"	****	<b>37,50</b>	<b>36,37</b>	<b>38,63</b>
1.1/4"	****	<b>31,50</b>	<b>30,55</b>	<b>32,45</b>
1"	****	<b>25,00</b>	<b>24,24</b>	<b>25,76</b>
3/4"	****	<b>19,00</b>	<b>18,42</b>	<b>19,58</b>
5/8"	****	<b>16,00</b>	<b>15,51</b>	<b>16,49</b>
1/2"	****	<b>12,50</b>	<b>12,11</b>	<b>12,89</b>
		<b>10,00</b>	<b>9,69</b>	<b>10,31</b>
3/8"	****	<b>9,50</b>	<b>9,20</b>	<b>9,80</b>
5/16"	****	<b>8,00</b>	<b>7,75</b>	<b>8,25</b>
1/4"	****	<b>6,30</b>	<b>6,10</b>	<b>6,50</b>
3,5	3,5	<b>5,60</b>	<b>5,42</b>	<b>5,78</b>
		<b>5,00</b>	<b>4,84</b>	<b>5,16</b>
4	4	<b>4,75</b>	<b>4,60</b>	<b>4,90</b>
5	5	<b>4,00</b>	<b>3,87</b>	<b>4,13</b>
6	6	<b>3,35</b>	<b>3,24</b>	<b>3,46</b>
7	7	<b>2,80</b>	<b>2,71</b>	<b>2,89</b>
8	8	<b>2,36</b>	<b>2,28</b>	<b>2,44</b>
10	9	<b>2,00</b>	<b>1,93</b>	<b>2,07</b>
12	10	<b>1,70</b>	<b>1,64</b>	<b>1,76</b>
14	12	<b>1,40</b>	<b>1,35</b>	<b>1,45</b>
16	14	<b>1,18</b>	<b>1,14</b>	<b>1,22</b>
18	16	<b>1,00</b>	<b>0,97</b>	<b>1,03</b>
*****	*****	<b>Abert. y</b>	<b>*****</b>	<b>*****</b>
20	20	<b>850</b>	<b>821</b>	<b>879</b>
25	24	<b>710</b>	<b>685</b>	<b>735</b>
30	28	<b>600</b>	<b>579</b>	<b>621</b>
35	32	<b>500</b>	<b>482</b>	<b>518</b>
40	35	<b>425</b>	<b>409</b>	<b>441</b>
45	42	<b>355</b>	<b>342</b>	<b>368</b>
50	48	<b>300</b>	<b>288</b>	<b>312</b>
60	60	<b>250</b>	<b>240,1</b>	<b>259,9</b>
70	65	<b>212</b>	<b>203,3</b>	<b>220,7</b>
80	80	<b>180</b>	<b>172,4</b>	<b>187,6</b>
100	100	<b>150</b>	<b>143,4</b>	<b>156,6</b>
120	115	<b>125</b>	<b>119,2</b>	<b>130,8</b>
140	150	<b>106</b>	<b>100,8</b>	<b>111,2</b>
170	170	<b>90</b>	<b>85,4</b>	<b>94,6</b>
200	200	<b>75</b>	<b>70,9</b>	<b>79,1</b>
230	250	<b>63</b>	<b>59,3</b>	<b>66,7</b>
270	270	<b>53</b>	<b>49,6</b>	<b>56,4</b>
325	325	<b>45</b>	<b>41,9</b>	<b>48,1</b>
400	400	<b>38</b>	<b>35,1</b>	<b>40,9</b>
500	500	<b>25</b>	<b>22</b>	<b>28</b>
635	635	<b>20</b>	<b>17</b>	<b>23</b>